OsteoMed® CFx™ Rigid Fixation System

Low Profile Orbital Floor Plate System

EN	Product Information and Instructions for Use	2
FR	Informations sur le produit et instructions d'utilisation	6
DE	Produktinformationen und Gebrauchsanweisung.....	11
IT	Informazioni sul prodotto e istruzioni per l'uso.....	16
PT	Informações sobre o produto e instruções de utilização	21
ES	Información sobre el producto e instrucciones de uso	26

Product Information and Instructions for Use

OsteoMed® CFX™ Rigid Fixation System Low Profile Orbital Floor Plate System

EN

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals in the European Union only.

Description

OSTEOMED CFX Rigid Fixation System is comprised of plates and screws used for fixation of craniofacial, maxillofacial, and mandibular fractures and reconstructions. The instruments are included to facilitate the placement of implants.

OSTEOMED CFX Rigid Fixation System can be used with the following:

- **OSTEOMED 2.0 Auto-Drive MMF Screws:** The OSTEOMED Auto-Drive MMF Screws are comprised of 2.0mm diameter screws used for temporary fixation to constrict fractured bone segments in the oral cavity. The system features screws in lengths from 8mm to 14mm.
- **OSTEOMED Low Profile Orbital Floor Plate System:** The OSTEOMED Low Profile Orbital Floor Plate System is offered with holes which can accept 1.2mm and 1.6mm standard and auto drive bone screws and 1.5mm and 1.9mm safety screws and are offered in a thickness ranging from .2mm to .5mm for use in the reconstruction of the floor and/or medial wall of the orbit.

Material

OSTEOMED CFX Rigid Fixation System: The screws are made from Titanium Alloy (ASTM F-136). The plates are made from Titanium Alloy (ASTM F-136) or commercially pure Titanium (ASTM F-67). The instrumentation is made from various grades of surgical stainless steel, anodized aluminum, and/or medical grade plastic.

OSTEOMED 2.0 Auto-Drive MMF Screws: The screws are made from Titanium Alloy (ASTM F-136). The wire is made from surgical grade stainless steel. The instrumentation is made from various grades of surgical grade stainless steel, anodized aluminum, and/or medical grade plastic.

OSTEOMED Low Profile Orbital Floor Plate System: The plates are made from commercially pure Titanium (ASTM F-67).

Intended Purpose

OSTEOMED CFX Rigid Fixation System is intended for fixation of craniomaxillofacial skeleton and bone fractures and reconstructions. Plates and screws are intended for single patient use only.

OSTEOMED 2.0 Auto-Drive MMF Screws are intended for temporary ligature and wire lock fixation for temporary constriction and stabilization of fractured bone segments in the oral cavity.

OSTEOMED Low Profile Orbital Floor Plate System is intended to be used in the reconstruction of the floor and/or medial wall of the orbit.

Indications for Use

Cranialmaxillofacial Indications Only

- **OSTEOMED CFX Rigid Fixation System** is indicated for fracture fixation during cranial and facial trauma reconstructions, orthognathic reconstruction, mandibular reconstruction, and surgery involving osteotomies and trauma of cranial or facial bones.
- **OSTEOMED Low Profile Orbital Floor Plate System** is indicated for the reconstructive treatment of orbital floor and/or medial wall trauma or bone excision in adult patients.

Maxillofacial Indications Only

- **OSTEOMED 2.0 Auto-Drive MMF Screws** are indicated for temporary ligature and wire lock fixation for temporary constriction and stabilization of fractured bone segments in the oral cavity in conjunction with primary fixation devices.

The intended clinical benefit of **OSTEOMED CFX Rigid Fixation System** is to provide fixation for patients with craniofacial, maxillofacial, and mandibular fractures as a result of trauma or genetic malformations of the craniomaxillofacial skeleton and bone. The indirect clinical benefit of the encompassed instrumentation is to facilitate implantation protocols for fracture fixation.

Contraindications

- Active or latent infection or marked inflammation of the treatment area
- Sepsis
- Insufficient quantity or quality of bone, osteoporosis, or in patients with certain metabolic diseases
- Patients with confirmed material sensitivity
- Patients who are unwilling or incapable of following post-operative care instructions
- 2.4mm Angled Locking Auto-Drive® screws are contraindicated in bi-cortical bone placement on the mandible

Warnings

- The treatment or implant may fail, including sudden failure, because of:
 - Loose fixation and/or loosening
 - Stress, including stress from inappropriate bending of the implant during surgery
 - Stress concentrations
 - Stress of weight bearing, load bearing, or excessive activity
 Failure is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. Failure is more likely if the patient does not follow post-operative care instructions.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant
- Instrument breakage or damage, as well as tissue damage, may occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use, or unintended use.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Fractured implants should be removed from patients during surgery. If unable to remove, notify the patient.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.

Cautions

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices in accordance with the Instructions for Use.
- Do not use sterile products past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Use of chemical disinfection may leave residues that adversely affect steam sterilization.

- Steam penetration and device sterilization may be negatively affected by labeling that block tray steam holes.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.
- When using the drill guide, do not apply a side load on the drill. This may result in friction, which may generate a thermal burn. Axial loading should always be used.
- When using the transbuccal approach, ensure cheek retractor is used to protect soft tissue.
- **OSTEOMED 2.0 Auto-Drive MMF Screws** and wire are intended for temporary fixation and are to be implanted for a maximum of four weeks.
- Manipulation of maxillofacial bones or tissues can potentially result in cardiac complications such as ectopic beats, atrioventricular block, bradycardia, syncope, vomiting, and asystole. These are known complications during maxillofacial surgery and are not specifically related to any device.

Potential Adverse Effects

- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Implant fracture due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at www.acumed.net/ifu.
- Injury to user.

Note: Any serious incident that has occurred in relation to the device should be reported to Acumed (through customercomplaints@acumed.net or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

Target Population

The system is intended to be used on skeletally mature, adult patients with sufficient quantity and quality of bone who require surgical intervention to achieve fracture fixation or reconstruction of cranial or facial bones.

MRI Status:

Many Acumed implants have been evaluated for safety in the MR environment and are MR conditional. Consult our publication "Acumed Implants in the MR Environment" at <http://www.acumed.net/ifu> for more information.

Device Lifetime:

Once installed, implants are expected to provide fixation and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon or patient. Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

Summary of Safety and Clinical Performance (SSCP)

The SSCP for the implants may be obtained from the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search term(s): 08456940BUDIO14NT

Maintaining Device Effectiveness

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation products and techniques.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- Acumed devices should be used in an operating room environment.
- All Acumed plates, screws, and instrumentation may be required for each surgery. Failure to use dedicated, unique Acumed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure that could cause subsequent patient injury. Failed devices may require re-operation and removal.
- Inspect the Acumed implants prior to use. Do not attempt a surgical procedure with faulty, damaged, or suspected instruments or implants. Alternate fixation methods should be available intraoperatively.
- Instrument lifetime is affected by usage, handling, and processing. Assess multiple-use instruments for fitness before and after each procedure to confirm that they are in proper operating condition. Indications that a reusable instrument should not continue to be used include but are not limited to deformation, dulled sharps, cracks, breakage, missing components, corrosion, excessive wear, and restriction or seizure of moving parts. Instruments which are faulty, damaged, or suspected should not be used. They should be replaced or sent to Acumed for disposition and repair.
- When placing more than one screw, ensure that subsequent screw placement does not interfere with other screws. Insert the second screw on the opposite side of the fracture or osteotomy site, and then all remaining screws, following the outlined procedures.
- Depth gauge marking tolerance: $\pm 0.25\text{mm}$.

Instructions for Use – OSTEOMED 2.0 Auto-Drive MMF Screws

1. Insert the screw in a TaperLock™ screwdriver and drive into the bone at a 90° angle using moderate pressure until the head is flush with the surface of the bone. Higher torque may be required to fully engage the threads than when using a normal screw with a pilot drill. Avoid multiple insertions of driver into the same screw to maintain self-retention feature of screw to driver.

Note: In high density bone, pilot drilling may be necessary. Drill Color Bands Denotes Size: 1.2mm drills are Yellow, 1.6mm drills are Blue and 2.4mm drills are Green. Speed and torque parameters must be in accordance with the Power System Instruction for Use. Use irrigation when pilot drilling.

Instructions for Use – Low Profile Orbital Floor Plate System

1. Expose the fracture site and retract soft tissue.
2. Select the orbital floor plate to be used based upon fracture location and patient anatomy.
3. If the orbital floor plate requires modification by cutting, plate cutter (P/N 220 0077) may be used to trim unneeded material.
4. Use deburring file (P/N 220 0584) if desired to remove any sharp edges.
5. If contouring is required prior to implantation, plate benders (P/Ns 220 0046 & 220 0049) may be used to pre bend the plate as needed.
Note: To prevent weakening of the plate, excessive and repeated bending should be avoided.
6. Insert and position plate with plate holding forceps (P/N 220 0027).
7. Fixate plate to bone using Acumed's patented Auto Drive™ self drilling screws, stop insertion of the screw once the screw is fully seated on the plate.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- Implants are single-use only and must not be reused. Reuse of single-use devices may result in transmission of infectious material from one patient to another which could result in serious injury or death.
- Contaminated implants must not be cleaned/reprocessed and reused. Any implant that has been contaminated by blood, tissue, and/or other bodily fluids/matter should never be used again and should be handled according to hospital protocol.
- Implants which are new, unused, and have never been contaminated by blood, tissue, and/or other bodily fluids/matter may be cleaned.
- Acumed recommends the following cleaning and sterilization instructions for implants and Instruments:

Manual Cleaning

IMPORTANT: Manual Cleaning steps are mandatory. Manual cleaning must be performed prior to performing automated cleaning.

1. Rinse the articles to be cleaned under running cool tap water to remove visible soil.
2. Use a syringe to flush water through any cracks, crevices, lumens, and hard to reach areas. During rinsing, actuate the articles to ensure thorough rinsing.
3. Prepare a neutral enzymatic cleaning solution per the manufacturer's instructions. Presoak soiled implants and instruments for a **minimum of 1 minute** in the enzymatic solution.

4. After soaking, thoroughly brush the articles beneath the surface of the cleaning solution using a soft bristled brush, paying close attention to all hard-to-reach areas until all evidence of soil is removed. Using a lumen brush or similar type brush, brush each lumen a minimum of 5 times, ensuring that the lumen brush is passed completely through or to the bottom of the entire lumen during brushing. Use twisting actions where possible. Actuate the articles while brushing to clean mated surfaces and movable parts.
5. Rinse implants and instruments under cool running water for at least 30 seconds to flush debris out of any tight spaces. Use a syringe to flush any lumen or mated surfaces.
6. Drain excess water from the article and dry using a clean, soft cloth and/or filtered pressurized air at ≤40 psi. Visually inspect each article without magnification for visible soil, deterioration, or loss of function. If soiled, repeat cleaning process. If the device is visibly deteriorated or is unable to function, return to Acumed.
7. Proceed to the automated cleaning steps -OR- Return the instruments to their sterilization case and steam sterilize according to Sterilization instructions.

Automated Cleaning

The following automated cleaning steps are optional. Automated cleaning using a washer/disinfector, according to the following steps, may be performed following the manual cleaning steps listed above.

8. Prepare a cleaning neutral enzymatic solution with the nominal detergent concentration per the washer/disinfector manufacturer's instructions.
9. Place implants within the appropriate Acumed organizer. Remove all trays from their cases prior to processing. Place the organizers and implants in a suitable washer/disinfector basket. Place instruments in a suitable washer/disinfector basket. Minimally invasive surgery (MIS) injector or irrigation ports may be affixed to cannulated instruments so cleaning spray forces water through the instrument lumens. Using the cleaning solution, process implants and instruments through a standard washer/disinfector cycle in compliance with EN ISO 15883-1 and EN ISO 15883-2 or equivalent national standards. Dry for the nominal time recommended by the washer/disinfector manufacturer. Do not use a lubrication phase.
10. Visually inspect each article under normal lighting and without magnification for residual water in blind holes, visible deterioration, or loss of function. Pay particular attention to all challenging areas. If the device is visibly deteriorated or is unable to function, return to Acumed.
11. Return the implants and instruments to their sterilization case and steam sterilize according to Sterilization instructions.

Thermal Disinfection

For thermal disinfection, remove trays from outer case and place separately in the washer/disinfector's rack during automated cleaning and thermal disinfection. Utilize a thermal disinfection program consistent with A0 value > 3000 or perform thermal disinfection at 90°C (maximum temperature 95°C) for at least 5 minutes.

Pre-Sterilization Inspection

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas. Reprocess instruments that are not clean and replace instruments that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace affected devices.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect cutting edges under magnification. Replace worn instruments e.g., dull, chipped, cracked, rolled, or otherwise deformed. Run a cotton cloth over the edge to help detect chipping and cracking.
- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.
- Lubricate (instrument milk) surgical instruments to increase useful life. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.

Sterility

- Product is supplied **NON-STERILE unless expressly labeled as STERILE**.
- Select plates and screws are available sterile packaged (Gamma Sterilized). **DO NOT USE IF STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and dry devices prior to sterilization per standard hospital procedures.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017, the Association for the Advancement of Medical Instrumentation's (AAMI's) "Comprehensive guide to steam sterilization and sterility assurance in health care facilities."
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of Acumed implant systems, the following parameters in the table below should be used. The table shows the minimum parameters validated to achieve a required Sterility Assurance Level (SAL) of 10⁻⁶ for the system.

If sterilizing the 210-3018 mini 4 hole straight plate, rigid, within the CFx Rigid Fixation System - Aluminum Tray, the plate will be placed in the locations marked for the 210-0030 plate.

	The CFx Rigid Fixation System- Aluminum Tray, is validated for steam sterilization using one of the below sets of parameters.		The CFx Rigid Fixation System- Plastic Tray, is validated for steam sterilization using one of the below sets of parameters.		The CFx Rigid Fixation System- Aluminum Tray, is validated for steam sterilization using one of the below sets of parameters.	
Pre-Vacuum Steam Sterilization	CFx Rigid Fixation System – Aluminum Tray*		CFx Rigid Fixation System – Plastic Tray		CFx Rigid Fixation System with Low Profile Orbital Floor Plate System – Aluminum Tray	
Temperature:	132°C	134°C	132°C	132°C	132°C	134°C
Full Cycle Time:	4 minutes	3 minutes	4 minutes	4 minutes	4 minutes	3 minutes
Configuration 1:	Wrapped tray		Wrapped tray		Wrapped tray	
Wrapped Dry Time Minimum	70 minutes	70 minutes	70 minutes	70 minutes	70 minutes	70 minutes
Wrapped Configuration Details:	Wrap in two layers of 1- ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques		Wrap in two layers of 1- ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques		Wrap in two layers of 1- ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques	
Configuration 2:	Rigid Containers	N/A	N/A		N/A	
Rigid Container Dry Time Minimum	30 minutes	N/A	N/A		N/A	
Rigid Sterilization Containers Configuration Details:	Remove modules and trays from CFx Rigid Fixation System Case. Load Full Length RSC (e.g., Aesculap JN442 base with JK489 lid) with two (2) plate modules on the bottom level and two (2) instrument trays stacked on top.	N/A	N/A		N/A	

*If sterilizing P/N 210-3018- mini 4-hole straight plate, rigid, within the CFx Rigid Fixation System- Aluminum Tray, the plate will be placed in the locations marked for the P/N 210-0030 plate.

Note: Biological indicator, *G. stearotherophilus*, was used in sterilization validation.

Post-Sterilization Inspection

- Do not store or use sterile devices if they are not dry.
- Inspect the sterile barrier for signs of damage. Do not use the product if the sterile barrier has been compromised.

Storage

After steam sterilization, sets should be stored under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017. Refer to sterilization wrap or rigid container manufacturer’s IFU for limits on sterile product storage time and storage requirements for temperature and humidity. Devices should be stored in an area that provides protection from dust, pests, and temperature/humidity extremes.

Sterile packaged implants should be stored at controlled room temperature out of direct sunlight and should be stored in an area that provides protection from dust and pests. Product package should be inspected prior to use for signs of damage or tampering.

Safe Disposal

Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling biohazardous material and disposal of sharp materials.

Symbols and Definitions					
Note: This is a general list of symbols. Not all symbols apply to all systems.					
Symbol	Description	ISO 15223-1	Symbol	Description	ISO 15223-1
	Manufacturer	5.1.1		Authorized Representative in the European Community	5.1.2
	Date of Manufacture	5.1.3		Use-by Date	5.1.4
	Batch Code (Lot Number)	5.1.5		Catalogue Number	5.1.6
	Serial number	5.1.7		Sterilized using irradiation	5.2.4
	Do Not Resterilize	5.2.6		Non-Sterile	5.2.7
	Do not use if sterile package is damaged	5.2.8		Single Sterile Barrier System	5.2.11
	Double Sterile Barrier System	5.2.12		Single Sterile Barrier System with Protective packaging inside	5.2.13
	Single Sterile Barrier System with Protective packaging inside	5.2.14		Single Use Only	5.4.2
	Medical Device	5.7.7	 <small>www.acumed.net/ifu</small>	Consult Instructions for Use	5.4.3
	MR Conditional (ASTM F2503)			CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment	
	Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician. (U.S. 21 CFR 801.109)				

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CUSTOMER SERVICE (OUTSIDE USA):

+1.888.627.9957

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-014204 Rev A | 07-2024 | © 2024 Acumed® LLC

OSTEOMED® is a registered trademark of Acumed, LLC.

Subsidiaries:

Acumed Ltd

Huebner House

The Fairground

Andover

Hampshire UK SP11 0QN

Tel: +44 1264 774450

Acumed GmbH

Fuhlsbuttler Stasse 300

22307 Hamburg

Deutschland

Tel: +49-40 947 82 093

Acumed Iberica C/

Procion 1 Edificio

Oficor 28023

Madrid-Spain Tel:

+34 913516357

EC REP

EMERGO EUROPE

WESTERVOORTSEDIJK 60

6827 AT ARNHEM

THE NETHERLANDS

CE
2797