## OsteoMed® 2.0 Orthognathic System

### Rigid Fixation System

EN	Product Information and Instructions for Use .....	2
FR	Informations sur le produit et instructions d'utilisation.....	7
DE	Produktinformationen und Gebrauchsanweisung ...	12
IT	Informazioni sul prodotto e istruzioni per l'uso .....	17
PT	Informações sobre o produto e instruções de utilização .....	21
ES	Información sobre el producto e instrucciones de uso.....	25

# Product Information and Instructions for Use

## OsteoMed® 2.0 Orthognathic System

### Rigid Fixation System

**These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals in the European Union only.**

#### **Description**

The **OSTEOMED 2.0 Orthognathic System** is comprised of plates, screws and instrumentation used for standard orthognathic maxillary and mandibular advancements.

#### **Intended Purpose**

**OSTEOMED 2.0 Orthognathic System** is intended for fixation of the craniomaxillofacial skeleton and bone fractures and reconstructions. Plates and screws are intended for single patient use only. Plates are intended for non-load bearing applications. The instruments are included to facilitate the placement of implants.

#### **Material**

The screws are made from Titanium Alloy (ASTM F-136). The plates are made from commercially pure Titanium (ASTM F-67). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade plastic. The chemical composition of Acumed implants is described in detail in the available Metal Sensitivity Statement document.

#### **Clinical Indications**

The **OSTEOMED 2.0 Orthognathic System** is indicated for fracture fixation in skeletally mature adults with mandibular fractures, maxillofacial fractures, or orthognathic needs as a result of trauma or genetic malformations of the craniomaxillofacial skeleton and bone.

The intended clinical benefit of the **OSTEOMED 2.0 Orthognathic System** is to provide fracture fixation for patients with mandibular fractures, maxillofacial fractures, or orthognathic needs as a result of trauma or genetic malformations of the craniomaxillofacial skeleton and bone. The indirect clinical benefit of the encompassed instrumentation is to facilitate implantation protocols for fracture fixation.

#### **Contraindications**

- Active or latent infection or marked inflammation of the treatment area
- Sepsis
- Insufficient quantity or quality of bone, osteoporosis, or in patients with certain metabolic diseases
- Patients with confirmed material sensitivity
- Patients who are unwilling or incapable of following post-operative care instructions

#### **Warnings**

- Treatment or implant may fail, including sudden failure, because of:
  - Loose fixation and/or loosening
  - Stress, including stress from inappropriate bending of the implant during surgery
  - Stress concentrations
  - Stress of weight bearing, load bearing, or excessive activityFailure is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. Failure is more likely if the patient does not follow post-operative care instructions.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Instrument breakage or damage, as well as tissue damage, may occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use, or unintended use.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Fractured implants should be removed from the patient during surgery. If unable to remove, notify the patient.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.

**Caution**

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices and instruments in accordance with the Instructions for Use.
- Do not use the sterile product past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Use of chemical disinfection may leave residues that adversely affect steam sterilization.
- Steam penetration and device sterilization may be negatively affected by labeling that block tray steam holes.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.
- Manipulation of maxillofacial bones or tissues can potentially result in cardiac complications such as but not limited to ectopic beats, atrioventricular block, bradycardia, syncope, vomiting, and asystole. These are known complications during maxillofacial surgery and are not specifically related to any device.

**Potential Adverse Effects**

- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Implant fracture due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur.
- Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document “Metal Sensitivity Statement” at [www.acumed.net/ifu](http://www.acumed.net/ifu).
- Injury to user.

**Target Population**

The system is intended to be used on skeletally mature, adult patients with sufficient quantity and quality of bone who require surgical intervention to achieve fracture fixation or reconstruction.

**MRI Status:**

Many Acumed implants have been evaluated for safety in the MR environment and are MR conditional. Consult our publication “Acumed Implants in the MR Environment” at [www.acumed.net/ifu](http://www.acumed.net/ifu) for more information.

**Device Lifetime:**

Once installed, implants are expected to provide fixation and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon or patient.

Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

**Summary of Safety and Clinical Performance (SSCP)**

The SSCP for the implants may be obtained from the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search term(s): 08456940BUDIO5LY.

Note: Any serious incident that has occurred in relation to the device should be reported to Acumed (through [customercomplaints@acumed.net](mailto:customercomplaints@acumed.net) or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

**Maintaining Device Effectiveness**

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation products and techniques.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- Acumed devices should be used in an operating room environment.
- All Acumed plates, screws, and instrumentation may be required for each surgery. Failure to use dedicated, unique Acumed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Inspect the Acumed implants prior to use. Do not attempt a surgical procedure with faulty, damaged, or suspected instruments or implants. Alternate fixation methods should be available intraoperatively.
- Instrument lifetime is affected by usage, handling, and processing. Assess multiple-use instruments for fitness before and after each procedure to confirm that they are in proper operating condition. Indications that a reusable instrument should not continue to be used include but are not limited to deformation, dulled sharp material, cracks, breakage, missing components, corrosion, excessive wear, and restriction or seizure of moving parts. Instruments which are faulty, damaged, or suspected should not be used. They should be replaced or sent to Acumed for disposition and repair.

**Instructions for Use, Zygomatic and Pyriform plate tabs**

1. Tabs are to be used for the identification of plates. They can also be used to aid in placement of plate on bone.
2. Tabs **must** be broken from plate once fixation is achieved. Bend tabs in opposite directions until broken from the plate.

**Instructions for Use, Depth Gauge**

1. If using bicortical fixation, slide the end of depth gauge through drilled pilot hole to measure screw length. Depth gauge marking tolerance:  $\pm 0.76\text{mm}$ .
2. Keeping depth gauge off the center of the pilot hole, pull back until the hook catches the lingual cortex.
3. Slide top sleeve down depth gauge until it rests on the plate or buccal cortex for positional screw placement.
4. Location of sleeve on measurement scale will indicate screw length.

**Cleaning**

- **Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.**

- Compliance is required with the equipment manufacturer’s user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- Implants are single-use only and must not be reused. Reuse of single-use devices may result in transmission of infectious material from one patient to another which could result in serious injury or death.
- Contaminated implants must not be cleaned/reprocessed and reused. Any implant that has been contaminated by blood, tissue, and/or other bodily fluids/matter should never be used again and should be handled according to hospital protocol.
- Implants which are new, unused, and have never been contaminated by blood, tissue, and/or other bodily fluids/matter may be cleaned.
- Acumed recommends the following cleaning and sterilization instructions for implants and Instruments:

**Manual Cleaning**

**IMPORTANT: Manual Cleaning steps are mandatory. Manual cleaning must be performed prior to performing automated cleaning.**

1. Rinse the articles to be cleaned under running cool tap water to remove visible soil.
2. Use a syringe to flush water through any cracks, crevices, lumens, and hard to reach areas. During rinsing, actuate the articles to ensure thorough rinsing.
3. Prepare a neutral enzymatic cleaning solution per the manufacturer’s instructions. Presoak soiled implants and instruments for a minimum of 1 minute in the enzymatic solution.
4. After soaking, thoroughly brush the articles beneath the surface of the cleaning solution using a soft bristled brush, paying close attention to all hard-to-reach areas until all evidence of soil is removed. Using a lumen brush or similar type brush, brush each lumen a minimum of 5 times, ensuring that the lumen brush is passed completely through or to the bottom of the entire lumen during brushing. Use a twisting action where possible. Actuate the articles while brushing to clean mated surfaces and movable parts.
5. Rinse implants and instruments under cool running water for at least 30 seconds to flush debris out of any tight spaces. Use a syringe to flush any lumen or mated surfaces.
6. Drain excess water from the article and dry using a clean, soft cloth and/or filtered pressurized air at ≤40 psi. Visually inspect each article without magnification for visible soil, deterioration, or loss of function. If soiled, repeat cleaning process. If the device is visibly deteriorated or is unable to function, return to Acumed.
7. Proceed to the automated cleaning steps -OR- Return the instruments to their sterilization case and steam sterilize according to Sterilization instructions.

**Automated Cleaning**

The following automated cleaning steps are optional. Automated cleaning using a washer/disinfector, according to the following steps, may be performed following the manual cleaning steps listed above.

8. Prepare a cleaning neutral enzymatic solution with the nominal detergent concentration per the washer/disinfector manufacturer’s instructions.
9. Place implants within the appropriate Acumed organizer. Remove all trays from their cases prior to processing. Place the organizers and implants in a suitable washer/disinfector basket. Place instruments in a suitable washer/disinfector basket. If possible, affix minimally invasive surgery (MIS) injector or irrigation ports to instruments so cleaning spray forces water and detergent through the instrument lumens. Using the cleaning solution, process implants and instruments through a standard washer/disinfector cycle in compliance with EN ISO 15883-1 and EN ISO 15883-2 or equivalent national standards. Dry for the nominal time recommended by the washer/disinfector manufacturer. Do not use a lubrication phase.
10. Visually inspect each article under normal lighting and without magnification for residual water in blind holes, visible deterioration, or loss of function. Pay particular attention to all challenging areas. If the device is visibly deteriorated or is unable to function, return to Acumed.
11. Return the implants and instruments to their sterilization case and steam sterilize according to Sterilization instructions.

**Thermal Disinfection:**

For thermal disinfection, remove trays from outer case and place separately in the washer/disinfector’s rack during automated cleaning and thermal disinfection. Cannulated instruments should be removed from trays and connected to an MIS injector unit or irrigation line to ensure adequate flow through the cannulation. Utilize a thermal disinfection program consistent with A0 value > 3000 or perform thermal disinfection at 90 °C (maximum temperature 95°C) for at least 5 minutes.

**Sterility**

All **OSTEOMED 2.0 Orthognathic System** plates, screws and instruments are packaged **NON-STERILE** and **MUST** be sterilized prior to surgical use.

<i>The 2.0 Orthognathic System is validated for steam sterilization using one of the below sets of parameters</i>		
<b>Pre-Vacuum Steam Sterilization</b>	<b>2.0 Orthognathic System</b>	
Temperature:	132°C	134°C
Full Cycle Time:	4 minutes	3 minutes
Dry Time:	55 minutes	55 minutes
Configuration:	Wrapped Tray	Wrapped Tray
Configuration Details:	Individually wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques.	Individually wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques.
Note: Biological indicator, <i>G. stearothermophilus</i> , was used in sterilization validation.		

Do not exceed 135°C, to avoid compromising functions of polymeric instrumentation.

**Storage**

Storage of Sets After Steam Sterilization:


Items should be stored under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017. Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.

Storage of Packaged Non-Sterile and Sterile Product:

Sterile devices should be stored in an area that provides protection from dust, pests, and temperature/ humidity extremes.

**Safe Disposal**

Devices may be contaminated with biohazardous materials and/or sharp. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling biohazardous material and disposal of sharp material.

<b>Symbols and Definitions</b>					
Note: This is a general list of symbols. Not all symbols apply to all systems.					
Symbol	Description	ISO 15223-1	Symbol	Description	ISO 15223-1
	Manufacturer	5.1.1		Authorized Representative in the European Community	5.1.2
	Date of Manufacture	5.1.3		Use-by Date	5.1.4
	Batch Code (Lot Number)	5.1.5		Catalogue Number	5.1.6
	Serial number	5.1.7		Sterilized using irradiation	5.2.4
	Do Not Resterilize	5.2.6		Non-Sterile	5.2.7
	Do not use if sterile package is damaged	5.2.8		Single Sterile Barrier System	5.2.11
	Double Sterile Barrier System	5.2.12		Single Sterile Barrier System with Protective packaging inside	5.2.13
	Single Sterile Barrier System with Protective packaging inside	5.2.14		Single Use Only	5.4.2
 <small>www.acumed.net/ifu</small>	Consult Instructions for Use	5.4.3		Medical Device	5.7.7
	MR Conditional (ASTM F2503)			CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment	
	Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician. (U.S. 21 CFR 801.109)				

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**ACUMED, LLC**  
**5885 NE CORNELIUS PASS ROAD**  
**HILLSBORO, OR 97124**

**CUSTOMER SERVICE (OUTSIDE USA):**  
**+1.888.627.9957**

*These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.*

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PKG1-014205-Rev A | Effective 07-2024 | © 2024 Acumed® LLC
Acumed® is a registered trademark of Acumed, LLC.
**Acumed Ltd**
**Acumed GmbH**