


Wrist Spanning Plates

The Acu-Loc 2 Wrist Spanning Plates are designed for complex distal radius fractures. It temporarily holds the wrist in distraction, providing ligamentotaxis to facilitate healing of the distal radius.



Key Features:

- ▶ Designed for second or third metacarpal placement
- ▶ Two plating options 171 mm and 188 mm lengths enables treatment of a wide range of patient anatomies
- ▶ 2.7 mm and 3.5 mm screw options (Ti or CoCr; Locking or Nonlocking)
- ▶ 2 mm of additional wrist distraction

*Than previous generation
†Data on file vs leading competitor

