



Acutrak® 3

Headless Compression Screw Solutions

Acutrak 3 offers **solutions for all 28 bones** in the foot, featuring **93 unique screws** available in **4 diameters**



Headless design allows subchondral implantation of the screws in and around articular regions with **minimal risk of impingement or soft tissue irritation**




**A Trusted Solution for Lower Extremities
Acutrak 3 Has Every Step Covered**

Acutrak 3

3+ Trusted Over
3 Million Times
Million

The Acutrak Headless Compression Screw, introduced in 1994, was the first fully threaded screw on the market with a continuously variable threaded pitch. Acutrak 3 is the latest innovation, featuring 93 unique screws, including a new 2.0 mm diameter Nano screw and extended lengths for the Micro, Mini, and Standard screws. These advancements expand the system's ability to address many indications historically treated with the Acutrak 2 System. Indicated for fractures, fusions, and osteotomies, the Acutrak family of screws is a versatile system designed to treat conditions in the forefoot, midfoot, and hindfoot.



Acutrak 3 Range of Sizes

Additional Sizes from Second Generations



4.0 Standard	3.5 Mini
2 mm increments 16-60 mm	2 mm increments 12-60 mm
2.5 Micro	2.0 Nano
1 mm increments 8-14 mm 2 mm increments 14-50 mm	1 mm increments 8-14 mm 2 mm increments 14-40 mm



3.5 Mini- Hallux Valgus



4.0 Standard – Tarsometatarsal Joint Fusion Involving Naviculocuneiform and second and third TMT joints




View More Acutrak 3 X-rays

go.acumed.net/AT3-Xrays

Expanded Solutions for Small Bones

Upgraded Nano and Micro lengths offer a broader range of sizes, delivering enhanced options for small bone procedures like hammer toe



Hallux Valgus Fixation, Redefined

With compression along its entire length, Acutrak 3 brings innovation and reliability to Hallux Valgus procedures

Contact Your Local Sales Rep for More Information



28 Operative Bones in the Foot & Ankle Acumed Can Treat Them All

Learn More: go.acumed.net/AT3



www.acumed.net

Oregon, USA Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
+1.888.627.9957

Texas, USA Campus
3885 Arapaho Road
Addison, TX 75001
+1.800.456.7779

Acumed Iberica Campus
C. Proci3n, 1
Edificio Oficor
28023 Madrid, Spain
+34.913.51.63.57

SPF10-29-A | Effective: 2025/02 | © 2025 Acumed® LLC

OsteoMed® LLC is a wholly owned subsidiary of Acumed LLC.
OsteoMed is a registered trademark of OsteoMed LLC.

Acumed® and Acutrak® are registered trademarks of Acumed LLC.